



Recombinant Mouse Acyl-protein thioesterase 2 (Lypla2)

Product Code	CSB-MP896227MO
Abbreviation	Lypla2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WTL7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MCGNTMSVPL LTDAATVSGA ERETAAVIFL HGLGDTGHSW ADALSTIRLP HVKYICPHAP RIPVTLNMKM VMPSWFDLMG LSPDAPEDEA GIKKAAENIK ALIEHEMKNQ IPANRIVLGG FSQGGALSLY TALTCPPLA GIVALSCWLP LHRNFPQAAAN GSAKDLAILQ CHGELDPMVP VRFGALTAEK LRTVVTPARV QFKTYPGVMH SSCPQEMAAV KEFLEKLLPP V
Source	Mammalian cell
Target Names	Lypla2
Protein Names	Recommended name: Acyl-protein thioesterase 2 Short name= APT-2 EC= 3.1.2.- Alternative name(s): Lysophospholipase 2 Lysophospholipase II Short name= LPL-II Short name= LysoPLA II Short name= mLyso II
Expression Region	1-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Lysophospholipases are enzymes that act on biological membranes to regulate the multifunctional lysophospholipids. There are alternatively spliced transcript variants described for this gene but the full length nature is not known yet.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.