



# Recombinant Arabidopsis thaliana Annexin D2 (ANN2)

<b>Product Code</b>	CSB-BP896377DOA
<b>Abbreviation</b>	ANN2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9XEE2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASLKVPSNV PLPEDDAEQL HKAFSGWGTN EKLIISILAH RNAAQRSLIR SVYAATYNED LLKALDKELS SDFERAVMLW TLDPPERDAY LAKESTKMFT KNNWVLVEIA CTRPALELIK VKQAYQARYK KSIEEDVAQH TSGDLRLLLL PLVSTFRYEG DDVNMMLARS EAKILHEKVS EKSYSDDDFI RILTTRSKAQ LGATLNHYNN EYGNAINKNL KEESDDNDYM KLLRAVITCL TYPEKHFEKV LRLSINKMGT DEWGLTRVVT TRTEVDMERI KEEYQRRNSI PLDRAIAKDT SGDYEDMLVA LLGHGDA
<b>Source</b>	Baculovirus
<b>Target Names</b>	ANN2
<b>Protein Names</b>	Recommended name: Annexin D2 Alternative name(s): AnnAt2
<b>Expression Region</b>	1-317
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.