



Recombinant Tursiops truncatus Interleukin-4 (IL4)

Product Code	CSB-BP896421TRF
Abbreviation	IL4
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9XS58
Product Type	Recombinant Protein
Immunogen Species	Tursiops truncatus (Atlantic bottle-nosed dolphin) (Delphinus truncatus)
Purity	≥85% (SDS-PAGE)
Sequence	HKCDVT LQEIIKTLNI LTAKKNLCME LPVEDVFATT KNTTEKETFC RAGTVLRHIY RHHKCFNQPL SGLHRNLSSM ANMTCSVNEA KKSTLKDFLE RLKMIMKEKY SKC
Source	Baculovirus
Target Names	IL4
Protein Names	Recommended name: Interleukin-4 Short name= IL-4 Alternative name(s): B-cell stimulatory factor 1 Short name= BSF-1 Lymphocyte stimulatory factor 1
Expression Region	25-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>