



Recombinant Probable 28S ribosomal protein S10, mitochondrial (Y37D8A.18)

Product Code	CSB-BP896445CXY
Abbreviation	Y37D8A.18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XWV5
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MLKLASSLRT GLISRSIRTL APTVNPAAEQQ QVQAVLPDKL YSSVEIEYRG HDKAVLKSYS SFLQQVCQHL EIPQGRLEVL PYIRWVQPAL RSKFVHKKYK LHYETRTHIS KLEILNVTGS TASTFLEYIQ RNIPPEGVGM R VGFTLQPLP LTIQKN
Source	Baculovirus
Target Names	mrps-10
Protein Names	Recommended name: Probable 28S ribosomal protein S10, mitochondrial Short name= MRP-S10 Short name= S10mt
Expression Region	1-156
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.