



Recombinant Drosophila melanogaster Mesencephalic astrocyte-derived neurotrophic factor homolog (Manf)

Product Code	CSB-YP896451DLU
Abbreviation	Manf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XZ63
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	LKEEDCEV CVKTVRRFAD SLDDSTKKDY KQIETAFKKF CKAQKNKEHR FCYYLGGLEE SATGILNELS KPLSWSMPAE KICEKLKKKD AQICDLRYEK QIDLNSVDLK KLVKRD LKKI LNDWDESCDG CLEKGFIDIKR IEELKPKYSR SEL
Source	Yeast
Target Names	Manf
Protein Names	Recommended name: Mesencephalic astrocyte-derived neurotrophic factor homolog Short name= DmMANF Alternative name(s): MANF/CDNF-like protein
Expression Region	23-173
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.