



Recombinant Human Serine/threonine-protein kinase 38-like (STK38L)

Product Code	CSB-EP896491HU-B
Abbreviation	STK38L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2H1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AMTAGTTTT FPMSNHTRER VTVAKLTLEN FYSNLILQHE ERETRQKKLE VAMEEEGLAD EEKCLRRSQH ARKETEFRLR KRTRLGLDDF ESLKVIGRGA FGEVRLVQKK DTGHIYAMKI LRKSDMLEKE QVAHIRAERD ILVEADGAWV VKMFYSFQDK RNLYLIMEFL PGGDMMTLLM KKDTLTEEET QFYISETVLA IDAIHQLGFI HRDIKPDNLL LDAKGHVKLS DFGLCTGLKK AHRTEFYRNL THNPPSDFSF QNMNSKRKAE TWKKNRRQLA YSTVGTDPDI APEVFMQTGY NKLCDWWSLG VIMYEMLIGY PFCSETPQE TYRKVMNWKE TLVFPPEVPI SEKAKDLILR FCIDSENRIK NSGVVEIKGH PFFEGVDWEH IRERPAAIPI EIKSIDDTSN FDDFPESDIL QPVPNTTEPD YKSKDWVFLN YTYKRFEGLT QRGSIPTYMK AGKL
Source	E.coli
Target Names	STK38L
Protein Names	Recommended name: Serine/threonine-protein kinase 38-like EC= 2.7.11.1 Alternative name(s): NDR2 protein kinase Nuclear Dbf2-related kinase 2
Expression Region	2-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.