



# Recombinant Human General transcription factor 3C polypeptide 5 (GTF3C5)

<b>Product Code</b>	CSB-BP896534HU
<b>Abbreviation</b>	GTF3C5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y5Q8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AAEAADLGL GAAVPVELRR ERRMVCVEYP GVVRDVAKML PTLGGEEGVS RIYADPTKRL ELYFRPKDPY CHPVCANRFS TSSLLLRIRK RTRRQKGVLG TEAHSEVTFD MEILGIISTI YKFQGMDFQ YLAVHTEAGG KHTSMYDKVL MLRPEKEAFF HQELPLYIPP PIFSRLDAPV DYFYRPETQH REGYNNPPIS GENLIGLSRA RRPNAIFVN FEDEEVPKQP LEAAQTWRR VCTNPVDRKV EEELRKLFDI RPIWSRNAVK ANISVHPDKL KVLLPFIAYY MITGPWRSLW IRFGYDPRKN PDAKIYQVLD FRIRCGMKHG YAPSDLPVKA KRSTYNYSLP ITVKKTSSQL VTMHDLKQGL GPSGTSGARK PASSKYKLD SVYIFREGAL PPYRQMFYQL CDLNVEELQK IHRNDGAEN SCTERDGWCL PKTSDELRTD MSLMIRQTIR SKRPALFSSS AKADGGKEQL TYESGEDEED EEEEEEEED FKPSDGSENE METEILDYV
<b>Source</b>	Baculovirus
<b>Target Names</b>	GTF3C5
<b>Protein Names</b>	Recommended name: General transcription factor 3C polypeptide 5 Alternative name(s): TF3C-epsilon Transcription factor IIIC 63 kDa subunit Short name= TFIIIC 63 kDa subunit Short name= TFIIIC63 Transcription factor IIIC subun
<b>Expression Region</b>	2-519
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.