



# Recombinant Human eIF5-mimic protein 1 (BZW2)

<b>Product Code</b>	CSB-MP896544HU
<b>Abbreviation</b>	BZW2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y6E2
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNKHQKPVLT GQRFKTRKRD EKEKFEPTVF RDTLVQQLNE AGDDLEAVAK FLDSTGSRLD YRRYADTLFD ILVAGSMLAP GGTRIDDGDK TKMTNHCVFS ANEDHETIRN YAQVFNKLIR RYKYLEKAFE DEMKKLLLFL KAFSETEQTK LAMLSGILLG NGTLPATILT SLFTDSLVE GIAASFAVKL FKAWMAEKDA NSVTSSLRKA NLDKRILLELF PVNRQSVDHF AKYFTDAGLK ELSDFLRVQQ SLGTRKELQK ELQERLSQEC PIKEVVLYVK EEMKRNDLPE TAVIGLLWTC IMNAVEWNKK EELVAEQALK HLKQYAPLLA VFSSQGQSEL ILLQKVQEYC YDNIHFMKAF QKIVVLFYKA DVLSEEAILK WYKEAHVAKG KSVFLDQMKK FVEWLQNAEE ESESEGEEN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	BZW2
<b>Protein Names</b>	Recommended name: Basic leucine zipper and W2 domain-containing protein 2
<b>Expression Region</b>	1-419
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.