



Recombinant Mungbean yellow mosaic virus Capsid protein (AR1, AV1)

Product Code	CSB-BP896663MNAR
Abbreviation	AR1, AV1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9YPS5
Product Type	Recombinant Protein
Immunogen Species	Mungbean yellow mosaic virus (strain Vigna) (MYMV)
Purity	>85% (SDS-PAGE)
Sequence	MPKRNYDTAF STPMSNVRRR LTFDTPLSLP ATAGSVPASA KRRRWNTNRPM WRKPRYYRLY RSPDVPRGCE GPCKVQSFEA KHDISHVGKV ICVTDVTRGM GITHRVGKRF CVKSIWVTGK IWMDENIKTK NHTNTVMFKL VRDRRPFGTP QDFGQVFNMY DNEPSTATVK NDLRDRYQVV RKFQATVTGG QYASKEQAIV SKFYRVNNYV VYNHQEAAKY ENHTENALLL YMACTHASNP VYATLKIRIY FYDSISN
Source	Baculovirus
Target Names	AR1
Protein Names	Recommended name: Capsid protein Alternative name(s): Coat protein Short name= CP
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.