



Recombinant Drosophila melanogaster Protein mab-21-like (CG4766)

Product Code	CSB-MP896681DLU
Abbreviation	CG4766
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y106
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLVPPDMMAA QTRMVYQMNR FCTERVQARM YKTATAIREI CKIVQDILKE VELQEPRFIS SLVECNRFVE GVEVISPNEF EIVLYLNQMG VFNFVDDGTL PGCAVLKLSLSD GRKRSMWLVW EFITASGYLS SRKIRARFQT LVAQACDKSV YRDMVKMVGDTSEVKLRIRE RFVVQITPAF KCSGIWPRSA AHWPVPHLPW PHPNIVAEVK TEGFDLLSKE SVIMQNKNNN AASMEGDAWV LSFFEAENRL LQGGCRRRCL SMLKTLRDRH LDLPGNPISA YHLKNLLLYE CEKHPRDFEW DEGCIADRIN GIFLQLISCL QYRRCPHYFL PALDMFKGKS PSALEQAQKQ VWRLTRELTT NANAFEKL
Source	Mammalian cell
Target Names	CG4766
Protein Names	Recommended name: Protein mab-21-like
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.