



Recombinant Human Protein AAR2 homolog (AAR2)

Product Code	CSB-EP896704HU-B
Abbreviation	C20orf4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y312
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAVQMDPEL AKRLFFEGAT VVILNMPKGT EFGIDYNSWE VGPKFRGVKM IPPGIHFLHY SSVDKANPKE VGPRMGFFLS LHQRGLTVLR WSTLREEVDL SPAPESEVEA MRANLQELDQ FLGPYPYATL KKWISLTNFI SEATVEKLQP ENRQICAFSD VLPVLSMKHT KDRVGQNLPR CGIECKSYQE GLARLPEMKP RAGTEIRFSE LPTQMFPEGA TPAEITKHSM DLSYALETVL NKQFPSSPQD VLGELQFAFV CFLLGNVYEA FEHWKRLNL LCRSEAMMK HHTLYINLIS ILYHQLGEIP ADFFVDIVSQ DNFLTSTLQV FFSSACSIIV DATLRKKA EK FQAHLTKKFR WDFAAEPEDC APVVVELPEG IEMG
Source	E.coli
Target Names	AAR2
Protein Names	Recommended name: Uncharacterized protein C20orf4
Expression Region	2-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.