



Recombinant Human PCTP-like protein (STARD10)

Product Code	CSB-EP896706HU
Abbreviation	STARD10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y365
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEKLAASTEP QGPRPVLGRE SVQVPDDQDF RSFRSECEAE VGWNLTYSRA GVSVWVQAVE MDRTLHKIKC RMECCDVPAE TLYDVLHDIE YRKKWDSNVI ETFDIARLTV NADVGYYSWR CPKPLKNRDV ITLRSWLPMG ADYIIMNYSV KHPKYPPRKD LVRAVSIQTG YLIQSTGPKS CVITYLAQVD PKGSLPKWVV NKSSQFLAPK AMKKMYKACL KYPEWKQKHL PHFKPWLHPE QSPLPSLALS ELSVQHADSL ENIDESAVAE SREERMGGAG GEGSDDDTSL T
Source	E.coli
Target Names	STARD10
Protein Names	Recommended name: PCTP-like protein Short name= PCTP-L Alternative name(s): Antigen NY-CO-28 START domain-containing protein 10 Short name= StARD10 Serologically defined colon cancer antigen 28 StAR-related lipid transfer
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.