



Recombinant Human Protein slowmo homolog 2 (SLMO2)

Product Code	CSB-YP896708HU
Abbreviation	SLMO2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y3B1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKIWTSEHVF DHPWETVTTA AMQKYPNPMN PSVVGVDVLD RHIDPSGKLH SHRLLSTEWG LPSIVKSLIG AARTKTYVQE HSVVDPVEKT MELKSTNISF TNMVSVDERL IYKPHPQDPE KTVLTQEAI TVKGVSLSSY LEGLMASTIS SNASKGREAM EWWIHKLNAE IEELTASARG TIRTPMAAAA FAEK
Source	Yeast
Target Names	PRELID3B
Protein Names	Recommended name: Protein slowmo homolog 2
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.