



# Recombinant Human Coiled-coil domain-containing protein 41 (CCDC41)

<b>Product Code</b>	CSB-BP896730HU
<b>Abbreviation</b>	CCDC41
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y592
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDTFPNNFPP GGDSGLTGSQ SEFQKMLIDE RLRCEHHKAN YQTLKAEHTR LQNEHVKLQN ELKHLFNEKQ TQQEKLQLLL EELRGELVEK TKDLEEMKLQ ILTPQKLELL RAQIQQELET PMRERFRNLD EEEVEKYRAVY NKLRYEHTFL KSEFEHQKEE YARILDEGKI KYESEIARLE EDKEELRNQL LNVDLTKDSK RVEQLAREKV YLCQKLGLE AEVAELKAEK ENSEAQVENA QRIQVRQLAE MQATVRSLEA EKQSANLRAE RLEKELQSSS EQNTFLINKL HKAEREINTL SSKVKELKHS NKLEITDIKL ETARAKSELE RERNKIQSEL DGLQSDNEIL KAAVEHHKVL LVEKDRELIR KVQAAKEEGY QKLVLQDEK LELENRLADL EKMKVEHDVW RQSEKDQYEE KLRASQMAEE ITRKELQSVR LKLQQQIVTI ENAEKEKNEN SDLKQQISSL QIQVTSLAQS ENDLLNSNQM LKEMVERLQK ECRNFRSQAE KAQLEAEKTL EEKQIQWLEE KHKLHERITD REEKYNQAKE KLQRAAIAQK KRKSLHENKL KRLQEKVEVL EAKKEELETE NQVLNRQNVP FEDYTRLQKR LKDIQRRHNE FRSLILVPMN PPTASINPVS FQSSAMVPSM ELPFPHPMQE EQHQRELSLL RKRLEELETT QRKQLEELGS SGE
<b>Source</b>	Baculovirus
<b>Target Names</b>	CEP83
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 41 Alternative name(s): Renal carcinoma antigen NY-REN-58
<b>Expression Region</b>	1-693
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.