



# Recombinant Human Ubiquitin carboxyl-terminal hydrolase isozyme L5 (UCHL5)

<b>Product Code</b>	CSB-EP896738HU-B
<b>Abbreviation</b>	UCHL5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y5K5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTGNAGEWCL MESDPGVFTE LIKGFGRGA QVEEWSLEP ENFEKLPVH GLIFLFWQP GEEPAGSVVQ DSRLDTIFFA KQVINNACAT QAIVSVLLNC THQDVHLGET LSEFKEFSQS FDAAMKGLAL SNSDVIRQVH NSFARQQMFE FDTKTSAKEE DAFHFVSYVP VNGRLYELDGL REGPIDLGA CNQDDWISAV RPVIEKRIQK YSEGEIRFNL MAIVSDRMI YEQKIAELQR QLAEEEPMDT DQGNMMLSAI QSEVAKNQML IEEVQKLKR YKIENIRRKH NYLPPFIMELL KTLAHQQLI PLVEKAKEKQ NAKKAQETK
<b>Source</b>	E.coli
<b>Target Names</b>	UCHL5
<b>Protein Names</b>	Recommended name: Ubiquitin carboxyl-terminal hydrolase isozyme L5 Short name= UCH-L5 EC= 3.4.19.12 Alternative name(s): Ubiquitin C-terminal hydrolase UCH37 Ubiquitin thioesterase L5
<b>Expression Region</b>	1-329
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.