



Recombinant Human Sorting nexin-8 (SNX8)

Product Code	CSB-BP896744HU
Abbreviation	SNX8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y5X2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTGRAMDPLP AAVGAAAEA EADEEADPPA SDLPTPQAIE PQAIVQQVPA PSRMQMPQGN PLLLSHTLQE LLARDTVQVE LIPEKKGLFL KHVEYEVSSQ RFKSSVYRRY NDFVVFQEML LHKFPYRMVP ALPPKRMLGA DREFIEARRR ALKRFVNLVA RHPLFSEDVV LKLFLSFGSGS DVQNKLKESA QCVGDEFLNC KLATRAKDFL PADIQAQFAI SRELIRNIYN SFHKLRDRAE RIASRAIDNA ADLLIFGKEL SAIGSDTTPL PSWAALNSST WGSCLKQALKG LSVEFALLAD KAAQQGKQEE NDVVEKLNLF LDLLQSYKDL CERHEKGV LH KHQRALHKYS LMKRQMMSAT AQNREPESVE QLESRIVEQE NAIQTMELRN YFSLYCLHQE TQLIHVYLPL TSHILRAFVN SSIQGHKEMS KVVNDLRPKL SCLFAGPHST LTPPCSPPED GLCPH
Source	Baculovirus
Target Names	SNX8
Protein Names	Recommended name: Sorting nexin-8
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.