



Recombinant *Schizosaccharomyces pombe* Uncharacterized J domain-containing protein C1778.01c (SPBC1778.01c, SPBC30D10.01)

Product Code	CSB-EP896772SXV
Abbreviation	SPBC1778.01c, SPBC30D10.01
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y7I8
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSAGDSIQLF PPTTKQAIEA SFVPHGKISP LIKRAVEPVG PLFLAHARRQ RHGRTFSEDE RLEVKNKVQE EVKEESEDEE EDPAMLRADP KEWKQQDHYA VLGLSKYRYK ADTEQIKKAH LKKVLKHHPD KKAASGNIND DSFFKCIQKA YEILSDPVR RQFDSVDENA DVEPPESTTK ETTFFELWTPV FESEARFSKK QPVPSLGTIE STRAEVDNFY NFWYNFDSWR SFEYLDKDIP DDGESRDNKR FQEKKNRSER QKNKARDNAR LRNLVDTALA SDPRIKLFKE QEKAACAARK WEREAGAREA AAAAQKKKEE EERRAAEEAA AKASAAAANK KAKEDKKKAQ KRDKKVVKNA LKDFNYFSAT DVPSAEHVDS VLKDVDVIMS KLGEGELGQL AADINAEKAA GAASVQAVFD KFAKMFIERG SMSSADVVF AQ
Source	E.coli
Target Names	zuo1
Protein Names	Recommended name: Uncharacterized J domain-containing protein C1778.01c
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.