



Recombinant Human PRELI domain-containing protein 1, mitochondrial (PRELID1)

Product Code	CSB-BP896867HU
Abbreviation	PRELID1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y255
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	TEDIV HREVTPDQKL LSRLLTKTN RMPRWAERLF PANVAHSVYV LEDSIVDPQN QTMTTFTWNI NHARLMVVEE RCVYCVNSDN SGWTEIRREA WVSSSLFGVS RAVQEFGLAR FKSNTKTMK GFEYILAKLQ GEAPSKTLVE TAKEAKEKAK ETALAATEKA KDLASKAATK KQQQQQQFV
Source	Baculovirus
Target Names	PRELID1
Protein Names	Recommended name: PRELI domain-containing protein 1, mitochondrial Alternative name(s): 25 kDa protein of relevant evolutionary and lymphoid interest Px19-like protein
Expression Region	36-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.