



# Recombinant Human Dexamethasone-induced Ras-related protein 1 (RASD1)

<b>Product Code</b>	CSB-BP896868HU
<b>Abbreviation</b>	RASD1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y272
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKLAAMIKKM CPSDSELSIP AKNCYRMVIL GSSKVGKTAI VSRFLTGRFE DAYTPTIEDF HRKFYSIRGE VYQLDILDTS GNHPFPAMRR LSILTGDVFI LVFSLDNRDS FEEVQLRQQLD IGTKSCLKN KTKENVDPVPL VICGNKGDRD FYREVDQREI EQLVGDDPQR CAYFEISAKK NSSLDQMFRA LFAMAKLPSE MSPDLHRKVS VQYCDVLHKK ALRNKLLRA GSGGGGGDPG DAFGIVAPFA RRPSVHSDLM YIREKASAGS QAKDKERC
<b>Source</b>	Baculovirus
<b>Target Names</b>	RASD1
<b>Protein Names</b>	Recommended name: Dexamethasone-induced Ras-related protein 1 Alternative name(s): Activator of G-protein signaling 1
<b>Expression Region</b>	1-278
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a Ras-related protein that is stimulated by dexamethasone. The exact function of this gene is unknown, but it may play a role in dexamethasone-induced alterations in cell morphology, growth and cell-extracellular matrix interactions. In addition, studies of a similar rat protein suggest that it functions as as a novel physiologic nitric oxide (NO) effector. The gene product belongs to the Ras superfamily of small GTPases.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.