



Recombinant Human Putative RNA-binding protein Luc7-like 2 (LUC7L2)

Product Code	CSB-EP896884HU
Abbreviation	LUC7L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y383
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSAQAQMRAM LDQLMGTSRD GDTTRQRIKF SDDRVCKSHL LNCCPHDVLS GTRMDLGECL KVHDLALRAD YEIASKEQDF FFELDAMDHL QSFIADCDRR TEVAKKRLAE TQEEISAEVA AKAERVHELN EEIGKLLAKV EQLGAEGNVE ESQKVMDEVE KARAKKREAE EVYRNSMPAS SFQQQKLRVC EVCSAYLGLH DNDRRLADHF GGKLHLGFIE IREKLEELKR VVAEKQEKRN QERLKRREER EREEREKLRR SRSHSKNPKR SRSREHRRHR SRSMSRERKR RTRSKSREKR HRHRSRSSSR SRSRSHQRSR HSSRDRSRER SKRRSSKERF RDQDLASCDR DRSSRDRSPR DRDRKDKKRS YESANGRSED RRSSEEREAG EI
Source	E.coli
Target Names	LUC7L2
Protein Names	Recommended name: Putative RNA-binding protein Luc7-like 2
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.