



Recombinant Human Coiled-coil domain-containing protein 9 (CCDC9)

Product Code	CSB-BP896890HU
Abbreviation	CCDC9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y3X0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAATLDLKSK EEKDAELDKR IEALRRKNEA LIRRYQEIEE DRKKAELEGV AVTAPRKGRS VEKENVAVES EKNLGPSRRS PGTPRPPGAS KGGRTPPQQG GRAGMGRASR SWEGSPGEQP RGGGAGGRGR RGRGRGSPHL SGAGDTSISD RKSKEWEERR RQNIKMNNEE MEKIAEYERN QREGVLEPNP VRNFLDDPRR RSGPLEESER DRREESRRHG RNWGGPDFER VRCGLEHERQ GRRAGLGSAG DMTLSMTGRE RSEYLRWKQE REKIDQERLQ RHRKPTGQWR REWDAEKTDG MFKDGPVPAH EPSHRYDDQA WARPPKPPTF GEFLSQHKA ASSRRRRKSS RPQAKAAPRA YSDHDDRWET KEGAASPAPE TPQPTSPETS PKETPMQPPE IPAPAHRPPE DEGEENEGEE DEEWEDISED EEEEEIEVEE GDEEPAQDH QAPEAAPTGI PCSEQAHGVP FSPEEPLLEP QAPGTPSSPF SPPSGHQPVSD DWGEEVELNS PRTHLAGAL SPGEAWPFES V
Source	Baculovirus
Target Names	CCDC9
Protein Names	Recommended name: Coiled-coil domain-containing protein 9
Expression Region	1-531
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.