



# Recombinant Human Coiled-coil domain-containing protein 9 (CCDC9)

<b>Product Code</b>	CSB-EP896890HU-B
<b>Abbreviation</b>	CCDC9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y3X0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAATLDLKSK EEKDAELDKR IEALRRKNEA LIRRYQEIEE DRKKAELEGV AVTAPRKGRS VEKENVAVES EKNLGPSRRS PGTPRPPGAS KGGRTPPQQG GRAGMGRASR SWEGSPGEQP RGGGAGGRGR RGRGRGSPHL SGAGDTSISD RKSKEWEERR RQNIKMNNEE MEKIAEYERN QREGVLEPNP VRNFLDDPRR RSGPLEESER DRREESRRHG RNWGGPDFER VRCGLEHERQ GRRAGLGSAG DMTLSMTGRE RSEYLRWKQE REKIDQERLQ RHRKPTGQWR REWDAEKTDG MFKDGPVPAH EPSHRYDDQA WARPPKPPTF GEFLSQHKAE ASSRRRRKSS RPQAKAAPRA YSDHDDRWET KEGAASPAPE TPQPTSPETS PKETPMQPPE IPAPAHRPPE DEGEENEGEE DEEWEDISED EEEEEIEVEE GDEEPAQDH QAPEAAPTGI PCSEQAHGVP FSPEEPLLEP QAPGTPSSPF SPPSGHQPVV DWGEEVELNS PRTHLAGAL SPGEAWPFES V
<b>Source</b>	E.coli
<b>Target Names</b>	CCDC9
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 9
<b>Expression Region</b>	1-531
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.