



Recombinant Human Leucine-rich repeat-containing protein 42 (LRRC42)

Product Code	CSB-MP896899HU
Abbreviation	LRRC42
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y546
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSYYLSSSENH LDPGPIYMRE NGQLHMVNLA LDGVRSSLQK PRPFRLFPKG FSVELCMNRE DDTARKEKTD HFIFTYTREG NLRYSKSLF SLVLGFISDN VDHIDSLIGF PEQIAEKLFS AAEARQKfte PGAGLRALQK FTEAYGSLVL CSLCLRNRYL VISEKLEEIK SFRELTCLDL SCCKLGDEHE LLEHLTNEAL SSVTQLHLKD NCLSDAGVRK MTAPVRVMKR GLENLTLDDL SCNPEITDAG IGYLFsFRKL NCLDISGTGL KDIKTVKHKL QTHIGLVHSK VPLKEFDHSN CKTEGWADQI VLQWERVTAE AVKPRETSEP RAAQRfYgk RSRAEAPLKC PLADTHMNSS EKLQFYKEKA PDCHGPVLKH EAISSQESKK SKKRPFESE TEQNNSSQPS KQKYVCLAVE DWDLNLSY
Source	Mammalian cell
Target Names	LRRC42
Protein Names	Recommended name: Leucine-rich repeat-containing protein 42
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.