



# Recombinant Human Leucine-rich repeat flightless-interacting protein 2 (LRRFIP2)

<b>Product Code</b>	CSB-MP896921HU
<b>Abbreviation</b>	LRRFIP2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y608
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGTPASGRKR TPVKDRFSAE DEALSNIARE AEARLAAKRA ARAEARDIRM RELERQQKEY SLHSFDRKWG QIQKWLEDSE RARYSHRSSH HRPYLGVEDA LSIRSVGSHR YDMFKDRSSR LSSLNHSYSH SHGMKKRSSD SHKDLLSGLY FDQRNYSSLR HSKPTSAYYT RQSSSLYSDP LATYKSDRAS PTANSGLLRS ASLASLYNGG LYNPYGPRTP SECSYYSSRI SSARSSPGFT NDDTASIVSS DRASRGRRES VVSAADYFSR SNRRGSVSE VDDISIPDLS SLDEKSDKQY AENYTRPSSR NSASATTPLS GNSSRRGSGD TSSLIDPDTLS LSELRDIYDL KDQIQDVEGR YMQGLKELKE SLSEVEEKYK KAMVSNAQLD NEKNNLIYQV DTLKDVIEEQ EEQMAEFYRE NEEKSKELER QKHMCSVLQH KMEELKEGLR QRDELIEEKQ RMQQKIDTMT KEVFDLQETL LWKDKKIGAL EKQKEYIACL RNERDMLREE LADLQETVKT GEKHGLVIIP DGTPNGDVSH EPVAGAITVV SQEAAQVLES AGEGPLDVRL RKLAGEKEEL LSQIRKLLKQ LEEERQKCSR NDGTVGDLAQ LQNGSDLQFI EMQRDANRQI SEYKFKLSKA EQDITTLEQS ISRLEGQVLR YKTAANAENAEK VEDELKAEKR KLQRELRTAL DKIEEMEMTN SHLAKRLEKM KANRTALLAQ Q
<b>Source</b>	Mammalian cell
<b>Target Names</b>	LRRFIP2
<b>Protein Names</b>	Recommended name: Leucine-rich repeat flightless-interacting protein 2 Short name= LRR FLII-interacting protein 2
<b>Expression Region</b>	1-721
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.