



Recombinant Human Leucine-rich repeat flightless-interacting protein 2 (LRRFIP2)

Product Code	CSB-MP896921HU
Abbreviation	LRRFIP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y608
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGTPASGRKR TPVKDRFSAE DEALSNIARE AEARLAAKRA ARAEARDIRM RELERQQKEY SLHSFDRKWG QIQKWLEDSE RARYSHRSSH HRPYLGVEDA LSIRSVGSHR YDMFKDRSSR LSSLNHSYSH SHGMKKRSSD SHKDLLSGLY FDQRNYSSLR HSKPTSAYYT RQSSSLYSDP LATYKSDRAS PTANSGLLRS ASLASLYNGG LYNPYGPRTP SECSYSSRI SSARSSPGFT NDDTASIVSS DRASRGRRES VVSAADYFSR SNRRGSVVSE VDDISIPDLS SLDEKSDKQY AENYTRPSSR NSASATTPLS GNSSRRGSGD TSSLIDPPTS LSELRDIYDL KDQIQDVEGR YMQLKELKE SLSEVEEKYK KAMVSNAQLD NEKNNLIYQV DTLKDVIEEQ EEQMAEFYRE NEEKSKELER QKHMCSVLQH KMEELKEGLR QRDELIEEKQ RMQQKIDTMT KEVFDLQETL LWKDKKIGAL EKQKEYIACL RNERDMLREE LADLQETVKT GEKHGLVIIP DGTPNGDVSH EPVAGAITVV SQEAAQVLES AGEGPLDVRL RKLAGEKEEL LSQIRKLLKQ LEEERQKCSR NDGTVGDLAG LQNGSDLQFI EMQRDANRQI SEYKFKLSKA EQDITTLEQS ISRLEGQVLR YKTAANAEEK VEDELKAEKR KLQRELRTAL DKIEEMEMTN SHLAKRLEKM KANRTALLAQ Q
Source	Mammalian cell
Target Names	LRRFIP2
Protein Names	Recommended name: Leucine-rich repeat flightless-interacting protein 2 Short name= LRR FLII-interacting protein 2
Expression Region	1-721
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.