



# Recombinant Human Protein FAM169A (FAM169A)

<b>Product Code</b>	CSB-YP896936HU
<b>Abbreviation</b>	FAM169A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y6X4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAFPVDMLEN CSHEELNSA EDYMSDLRCG DPENPECFSL LNITIPISLS NVGFVPLYGG DQTQKILALF APEDSLTAVA LYLADQWWAI DDIVKTSVPS REGLKQVSTL GERVVLYVLN RIIYRKQEME RNEIPFLCHS STDYAKILWK KGEAIGFYSV KPTGSICASF LTQSYQLPVL DTMFLRKKYR GKDFGLHMLE DFVDSFTEDA LGLRYPLSSL MYTACKQYFE KYPGDHELLW EVEGVGHWYQ RIPVTRALQR EALKILALSQ NEPKRPMSE YGPASVPEYE ARTEDNQSSE MQLTIDSLKD AFASTSEGHD KTSVSTHTRS GNLRPKKIGK RFQDSEFSSS QGEDEKTSQT SLTASINKLE STARPSESSE EFLEEEPEQR GIEFEDESSD RDARPALETQ PQQEKQDGEK ESELEPMNGE IMDDSLKTSL ITEEDSTSE VLDEELKLQP FNSSEDSTNL VPLVVESSKP PEVDAPDKTP RIPDSEMLMD EGTSDEKGHM EEKLSLLPRK KAHLGSSDNV ATMSNEERSD GGFPNSVIAE FSEEPVSENL SPNTTSSLED QGEEGVSEPPQ ETSTALPQSS LIEVELEDVP FSQNAGQKNQ SEEQSEASSE QLDQFTQSAE KAVDSSSEI EVEVPVDRR NLRRKAKGKH GPAKKKAKLT
<b>Source</b>	Yeast
<b>Target Names</b>	FAM169A
<b>Protein Names</b>	Recommended name: Protein FAM169A
<b>Expression Region</b>	1-670
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.