



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C685.03 (SPBC685.03)

Product Code	CSB-EP896951SXV-B
Abbreviation	SPBC685.03
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y7L5
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	LPSPGIPSND KDLAIYPPVD ANNTTSQDVV IDVVLFGNGTE WISSRERFSL VSPGWASIYG PSGQSSVGSV LYSPILDYID NAQLLVLDDV SFGVFSSLAN SVTGYASNYT QSSSNSTSTM NSTGVSVSGGS VYPTNSTTNS SISWNSSTSA ATNTSSSSSS SSQSSVSVN SEIFSIFGLS QQYVNYSTSR LCVVGTPRAN MSTVSVTNG SAVSNYTVNT NGWTSSNFKC VDDVANIFG LDFYTAAVLS EVSILRSFAL CNATTSSSLF RQIASYGVIY SFHFSSSESG SFANLIGTNN YFMTDVKSSS VVIVQSETSC SINSASMSSN TTYFYWNSTS SLSSSVFTNT TSSSNSTNSS IPTTYPNST TYQNITTSYP WSQPVVNITD YLSDNGDGHF VLAGDGNQTI GDFYVMNWTT IASGEYLVPF NY
Source	E.coli
Target Names	SPBC685.03
Protein Names	Recommended name: Uncharacterized protein C685.03
Expression Region	21-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.