



# Recombinant Schizosaccharomyces pombe Pre-mRNA-splicing factor cwc25 (cwf25)

<b>Product Code</b>	CSB-EP896968SXV-B
<b>Abbreviation</b>	cwf25
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y805
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGGGDLNMKK SWHPLLMRNQ EKVWKDEQAH KEEMKRVEQL RREIEEERQL LELHRLQEAA GGKKRKDRVE WMYAVPNTNG PNRDSSEMEE YLLGRRRLDD LLKDKIEDQN NSLEKTEFIA LQNANSLQDT QAKLRLDPLL AIKQQEQKQL QTLMEKRKYS LDSDRKSKER RHRDRHHRSN QDRSRERSDN EQHSSDKREH SRRSYRNDRN NWRERTHNR YRHRDKYDSG YFKKHYYDDM RFDQGHFQDE RDLKKYVRTS RQYSRSPSPD FRTRNHQFHS RDSQPITQRH TDIESRLQKM QDNAKELDES RRKKIELLEK KERDEEQFLE KERRDTARKW DNQGDFIRNM RKEIYSGDSV SLADRVNSSR HNMLRP
<b>Source</b>	E.coli
<b>Target Names</b>	cwf25
<b>Protein Names</b>	Recommended name: Pre-mRNA-splicing factor cwc25 Alternative name(s): Cell cycle control protein cwf25 Complexed with cdc5 protein 25
<b>Expression Region</b>	1-376
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.