



Recombinant Probable alpha-galactosidase B (aglB)

Product Code	CSB-YP896972DOZ
Abbreviation	aglB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y865
Product Type	Recombinant Protein
Immunogen Species	Aspergillus niger
Purity	≥85% (SDS-PAGE)
Sequence	LVRP DGVGLTPALG WNSWNAYSCD IDADKIVTAA NEVVNLGLKD LGYEYINIDD CWSVKSGRNT TTKRIIPDPD KFPNGISGVA DQVHALGLKL GIYSSAGLTT CAGYPASLGY EEIDAQSFAE WGIDYLKYDN CGVPTNLTDQ YTYCVPDSTD GSNYPNGTCV NLDAAPQGY DWATSTTAKR YQRMRDALLS VNRTILYSLC DWGQADVNAW GNATGNSWRM SGDITATWSR IAEIANENSF LMNYANFWGY PDPDMLEVGN GNLTLPENRA HFALWAMMKA PLIIGTPLDS IDTSHLTILS NKPLLTFHQD AVIGRPAYPY KWGYNPDWTF DPEHPAEYWS GPTSSGEV FV LMLNSEGEVK TRSAVWEEVP ELKDRGTKKN SKEKKGFKVT DAWTGKDLGC VKDKYEVKLQ AHDVAVLVVG GQC
Source	Yeast
Target Names	aglB
Protein Names	Recommended name: Probable alpha-galactosidase B EC= 3.2.1.22 Alternative name(s): Melibiase B
Expression Region	17-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.