



Recombinant Aeropyrum pernix 50S ribosomal protein L22P (rpl22p)

Product Code	CSB-BP897021AYB
Abbreviation	rpl22p
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9YF76
Product Type	Recombinant Protein
Immunogen Species	Aeropyrum pernix (strain ATCC 700893 / DSM 11879 / JCM 9820 / NBRC 100138 / K1)
Purity	>85% (SDS-PAGE)
Sequence	MPRWGYSVKL RDESEVAKAV LRNVPVHPKI MAEVARAISG MRVDEARRYL RAVIEKREAV PFRRAHGKQA HRRGLADKWG WPVGRYPVKA ARYMLKLLDN VEANAANKNL DVERLKIIVH AAHKGITLKR WMPRAWGRAT PRNRVHSHIE IMVREV
Source	Baculovirus
Target Names	rpl22
Protein Names	Recommended name: 50S ribosomal protein L22P
Expression Region	1-156
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.