



Recombinant Mouse Cell division cycle 7-related protein kinase (Cdc7)

Product Code	CSB-EP897068MO-B
Abbreviation	Cdc7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Z0H0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MEEPMAFSSL RGS DRCPADD SLKKYEQSVK LSGIKRDIEE LCEAVPQLVN VFKIKDKIGE GTFSSVYLAT AQLQEGHEEK IALKHLIPTS HPMRIA AELQ CLTVAGGQDN VMGLKYCFRK NDHVVIAMPY LEHESFLDIL NSLSFQEVRE YMYNLFVALK RIHQFGIVHR DVKPSNFLYN RRLKKYALVD FGLAQGTRDT KIELLK FVQS EAQQEDCSR N KYHG VVG H KG LLSRPAPKTV DQQCTPKTSV KRSY TQVHIK QGKDGKERSV GLSVQRSVFG ERNFNIHSSI SHESPAEKLI KQSKTVDIIS RKLATK KTAI STKAMNSVMR ETARSCPAVL TCDCYGS DRV CSVCLSR RQQ VAPRAGTPGF RAPEVLTKCP DQTTAIDMWS AGVIFLSLLS GRYPFYKASD DLTALAQIMT IRGSRETIQA AKAFGKSVLC SKEVPAQDLR ALCERLRGLD STTPRSASGP PGNASYDPAA SKNTD HKASR VQAAQAQHSE DSLYKRDNDG YW SHPKDCTS NSEGWDSVPD EAYDLLDKLL DLNPASRITA EAALLHAFFK DMCS
Source	E.coli
Target Names	Cdc7
Protein Names	Recommended name: Cell division cycle 7-related protein kinase Short name= CDC7-related kinase Short name= muCdc7 EC= 2.7.11.1
Expression Region	1-564
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a cell division cycle protein with kinase activity that is critical for the G1/S transition. The yeast homolog is also essential for initiation of DNA replication as cell division occurs. Overexpression of this gene product may be associated with neoplastic transformation for some tumors. Multiple alternatively spliced transcript variants that encode the same protein have been detected.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.