



Recombinant *Drosophila melanogaster* NuA4 complex subunit EAF3 homolog (MRG15)

Product Code	CSB-MP897072DLU
Abbreviation	MRG15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y011
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGEVKPAKVE NYSTGTDANT LFVDGERVLC FHGPLIYEAK VLKTKPDATP VEYYIHYAGW SKNWDEWVPE NRVLKYNDN VKRRQELARQ CGERSKKDNK KGSAAKME QMRNESRAST PSKDSNTSQS TASSTPTTSA GPGSKSEAGS TGTTTTNSTA NSTTSRAHRK STQSTPSTAR PGTPSDKKED PAAAEETEEE GPVAPKKKRM SEQRPSLTGS DVAEKPLPPT TTPSTPTTEP APCVESEEAY AAKVEVKIKI PDELKHYLTD DWYAVVREHK LLELPAKVTV QQISEQYLAH KKSVKSTSAS KEVAINDVLD GIVEYFNVML GSQLLYKFER TQYADVMQKH PDTPLSELYG SFHLLRLFVR LGSMLSYSAL DQQSMQNLIT HVQDFLKFLV KNSSIFFSMS NFINVDPEYV RNAQ
Source	Mammalian cell
Target Names	MRG15
Protein Names	Recommended name: NuA4 complex subunit EAF3 homolog Alternative name(s): Protein MRG15
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.