



Recombinant Human Trafficking protein particle complex subunit 4 (TRAPPC4)

Product Code	CSB-YP897085HU
Abbreviation	TRAPPC4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y296
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAIFSVYVVN KAGGLIYQLD SYAPRAEAEK TFSYPLDLLL KLHDERVLVA FGQRDGIRVG HAVLAINGMD VNGRYTADGK EVLEYLGNPA NYPVSIRFGR PRLTSNEKLM LASMFHSLFA IGSQLSPEQG SSGIEMLETD TFKLHCYQTL TGIKFVVLAD PRQAGIDSLL RKIYEIYSDF ALKNPFYSLE MPIRCELFQD NLKLALAEVAE KAGTFGPGS
Source	Yeast
Target Names	TRAPPC4
Protein Names	Recommended name: Trafficking protein particle complex subunit 4 Alternative name(s): Hematopoietic stem/progenitor cell protein 172 Synbindin TRS23 homolog
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.