



# Recombinant Human Tudor and KH domain-containing protein (TDRKH)

<b>Product Code</b>	CSB-EP897096HU-B
<b>Abbreviation</b>	TDRKH
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y2W6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSTERTSWTS LSTIQKIALG LGIPASATVA YILYRRYRES REERLTFVGE DDIEIEMRVP QEAVKLIIGR QGANIKQLRK QTGARIDVDT EDVGDERVLL ISGFPVQVCK AKAAIHQILT ENTPVSEQLS VPQRSVGRRI GRGGETIRSI CKASGAKITC DKESEGTL L SRLIKISGTQ KEVAAAKHLI LEKVSEDEEL RKRIAHSAET RVPRKQPISV RREDMTEPGG AGEPALWKNT SSSMEPTAPL VTPPPKGGGD MAVVVSKEGS WEKPSDDSFQ KSEAQAIPEM PMFEIPSPDF SFHADEYLEV YVSASEHPNH FWIQIVGSR S LQLDKLVNEM TQHYENSVPE DLTVHVGDIV AAPLPTNGSW YRARVLGTLE NGNLDLYFVD FGDNGDCPLK DLRALRSDFL SLPFQAIECS LARIAPSGDQ WEEEALDEFD RLTHCADWKP LVAKISSYVQ TGISTWPKIY LYDTSNGKKL DIGLELVHKG YAIELPEDIE ENRAVPDMLK DMATETDASL STLLTETKKS SGEITHTLSC LSLSEAASMS GDDNLEDDYL L
<b>Source</b>	E.coli
<b>Target Names</b>	TDRKH
<b>Protein Names</b>	Recommended name: Tudor and KH domain-containing protein Alternative name(s): Tudor domain-containing protein 2
<b>Expression Region</b>	1-561
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.