



Recombinant Human Tudor and KH domain-containing protein (TDRKH)

Product Code	CSB-EP897096HU
Abbreviation	TDRKH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2W6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSTERTSWTS LSTIQKIALG LGIPASATVA YILYRRYRES REERLTFVGE DDIEIEMRVP QEAVKLIIGR QGANIKQLRK QTGARIDVDT EDVGDERVLL ISGFPVQVCK AKAAIHQILT ENTPVSEQLS VPQRSVGRRII GRGGETIRSI CKASGAKITC DKESEGTHLL SRLIKISGTQ KEVAAAKHLI LEKVSEDEEL RKRIAHS AET RVPRKQPISV RREDMTEPGG AGEPALWKNT SSSMEPTAPL VTPPPKGGGD MAVVVSKEGS WEKPSDDSFQ KSEAQAIPEM PMFEIPSPDF SFHADEYLEV YVSASEHPNH FWIQIVGSR S LQLDKLVNEM TQHYENSVPE DLTVHVGDIV AAPLPTNGSW YRARVLGTLE NGNLDLYFVD FGDNGDCPLK DLRALRSDFL SLPFQAIECS LARIAPSGDQ WEEEALDEFD RLTHCADWKP LVAKISSYVQ TGISTWPKIY LYDTSNGKKL DIGLELVHKG YAIELPEDIE ENRAVPDMLK DMATETDASL STLLTETKKS SGEITHTLSC LSLSEAASMS GDDNLEDDYL L
Source	E.coli
Target Names	TDRKH
Protein Names	Recommended name: Tudor and KH domain-containing protein Alternative name(s): Tudor domain-containing protein 2
Expression Region	1-561
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.