



Recombinant Human Ribosomal biogenesis protein LAS1L (LAS1L)

Product Code	CSB-EP897285HU-B
Abbreviation	LAS1L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y4W2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSWESGAGPG LGSQGMDLVW SAWYGKCVKG KGSLPLSAHG IVVAWLSRAE WDQVTVYLFC DDHKLQRYAL NRITVWRSRS GNELPLAVAS TADLIRCKLL DVTGGLGTDE LRLLYGMALV RFVNLISERK TKFAKVPLKC LAQEVNIPDW IVDLRHELTH KKMPHINDCR RGCYFVLDWL QKTYWCRQLE NSLRETWELE EFREGIEEED QEEDKNIVVD DITEQKPEPQ DDGKSTESDV KADGDSKGSE EVDSHCKKAL SHKELYERAR ELLVSYEEEQ FTVLEKFRYL PKAIKAWNNP SPRVECVLAE LKGVTCENRE AVLDAFLDDG FLVPTFEQLA ALQIEYEDGQ TEVQRGEGTD PKSHKNVDLN DVLVPKPFSSQ FWQPLLRLGLH SQNFTQALLE RMLSELPALG ISGIRPTYIL RWTVELIVAN TKTGRNARRF SAGQWEARRG WRLFNCSASL DWPRMVESCL GSPCWASPQL LRIIFKAMGQ GLPDEEQEKL LRICSIYTS GENS LVQEGS EASPIGKSPY TLDSLYWSVK PASSSFGSEA KAQQQEEQGS VNDVKEEKE EKEVLPDQVE EEEENDDQEE EEEDEDEDDED EEDRMEVGP FSTGQESPTA ENARLLAQKR GALQGS AWQV SSEDVRWDTF PLGRMPGQTE DPAELMLENY DTMYLLDQPV LEQRLEPSTC KTDTLGLSCG VGSGNCSNSS SSNFEGLLWS QGQLHGLKTG LQLF
Source	E.coli
Target Names	LAS1L
Protein Names	Recommended name: Ribosomal biogenesis protein LAS1L Alternative name(s): Protein LAS1 homolog
Expression Region	1-734
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.