



# Recombinant Human Ribosomal biogenesis protein LAS1L (LAS1L)

<b>Product Code</b>	CSB-EP897285HU-B
<b>Abbreviation</b>	LAS1L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y4W2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSWESGAGPG LGSQGMDLVW SAWYGKCVKG KGSLPLSAHG IVVAWLSRAE WDQVTVYLFC DDHKLQRYAL NRITVWRSRS GNELPLAVAS TADLIRCKLL DVTGGLGTDE LRLLYGMALV RFVNLISERK TKFAKVPLKC LAQEVNIPDW IVDLRHELTH KKMPHINDCR RGCYFVLDWL QKTYWCRQLE NSLRETWELE EFREGIEEED QEEDKNIVVD DITEQKPEPQ DDGKSTESDV KADGDSKGSE EVDSHCKKAL SHKELYERAR ELLVSYEEEQ FTVLEKFRYL PKAIKAWNNP SPRVECVLAE LKGVTCENRE AVLDAFLDDG FLVPTFEQLA ALQIEYEDGQ TEVQRGEGTD PKSHKNVDLN DVLVPKPFSQ FWQPLLRLGLH SQNFTQALLE RMLSELPALG ISGIRPTYIL RWTVELIVAN TKTGRNARRF SAGQWEARRG WRLFNCASL DWPRMVESCL GSPCWASPQL LRIIFKAMGQ GLPDEEQEKL LRICSIYTS GENSLVQEGS EASPIGKSPY TLDSLYWSVK PASSSFGSEA KAQQQEEQGS VNDVKEEKE EKEVLPDQVE EEEENDDQEE EEEDEDEDDED EEDRMEVGP FSTGQESPTA ENARLLAQKR GALQGS AWQV SSEDVRWDTF PLGRMPGQTE DPAELMLENY DTMYLLDQPV LEQRLEPSTC KTDTLGLSCG VGSGNCSNSS SSNFEGLLWS QGQLHGLKTG LQLF
<b>Source</b>	E.coli
<b>Target Names</b>	LAS1L
<b>Protein Names</b>	Recommended name: Ribosomal biogenesis protein LAS1L Alternative name(s): Protein LAS1 homolog
<b>Expression Region</b>	1-734
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.