



Recombinant Human AP-4 complex subunit sigma-1 (AP4S1)

Product Code	CSB-EP897288HU
Abbreviation	AP4S1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y587
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MIKFFLMV NK QGQTRL SKYY EHVDINKRTL LETEVIKSCL SRSNEQCSFI EYKDFKLIYR QYAALFIVVG VNDTENEMAI YEFIHNFVEV LDEYFSRVSE LDIMFNLDKV HIILDEMVLN GCIVETNRAR ILAPLLILDK MSES
Source	E.coli
Target Names	AP4S1
Protein Names	Recommended name: AP-4 complex subunit sigma-1 Alternative name(s): AP-4 adapter complex subunit sigma-1 Adapter-related protein complex 4 subunit sigma-1 Sigma-1 subunit of AP-4 Sigma-4-adaptin Short name= Sigma4-adaptin
Expression Region	1-144
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.