



# Recombinant Mouse Tyrosine-protein kinase TXK (Txk)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP025362MO-B  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P42682  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Mus musculus (Mouse)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MILSSYSSFQ SVLCCCCRC SVQKRQVRTQ ISLSREEELS EKHSQRQRPW<br>FAKLMGKTQS NRGGVQPSKR KPLPPLPQEP PDERIQVKAL YDFLPREPGN<br>LALKRAEEYL ILERCDPHWW KARDRFGNEG LIPSNYVTEN RLANLEIYEW<br>YHKNITRNQT ERLLRQEAKE GAFIVRDSRH LGSYTISVFT RARRHTQSSI<br>KHYQIKKND S GQWYITERHL FPSVPELIQY HQYNAAGLIS RLRYPILLG<br>SCLPATSGFS YEKWEIDPSE LAFVKEIGSG QFGVVHLGEW RAHIPVAIKA<br>INEGSMSEED FIEEAKVMMK LSHSRLVQLY GVCIQKPLY IVTEFMENG<br>LLDYLRRERKG QLQKALLLSM CQDICEGMAY LERSCYIHRD LAARNCLVSS<br>ACVVKISDFG MARYVLDDEY ISSSGAKFPV KWCPPEVFHF NKYSSKSDVW<br>SFGVLMWEVF TEGKMPFENK SNLQVVEAIS QGFRLYRPHL APMTIYRVMY<br>SCWHESPKGR PTFAELLQVL TEIAETW |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | Txk   |
| <b>Protein Names</b>     | Recommended name: Tyrosine-protein kinase TXK EC= 2.7.10.2<br>Alternative name(s): PTK-RL-18 Resting lymphocyte kinase  |
| <b>Expression Region</b> | 1-527   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |