



Recombinant Human Cytoplasmic dynein 1 light intermediate chain 1 (DYNC1LI1)

Product Code	CSB-MP897312HU
Abbreviation	DYNC1LI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y6G9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAVGRVGSF GSSPPGLSST YTGGPLGNEI ASGNGGAAAG DDEDGQNLWS CILSEVSTRS RSKLPAGKNV LLLGEDGAGK TSLIRKIQGI EEYKKGRGLE YLYLNVHDED RDDQTRCNVW ILGDLYHKG LLKFSLDAVS LKDTLVMLVV DMSKPWTALD SLQKWASVVR EHVDKLKIPP EEMKQMEQKL IRDFQEYVEP GEDFPASPQR RNTASQEDKD DSVVLPLGAD TLTHNLGIPV LVVCTKCDAI SVLEKEHDYR DEHFDFIQSH IRKFCLQYGA ALIYTSVKEN KNIDLVIKYI VQKLYGFPYK IPAVVVEKDA VFIPAGWDND KKIGILHENF QTLKAEDNFE DIITKPPVRK FVHEKEIMAE DDQVFLMKLQ SLLAKQPPTA AGRPVDASPR VPGGSPRTPN RSVSSNVASV SPIPAGSKKI DPNMKAGATS EGLANFFNS LLSKKTGSPG GPGVSGGSPA GGAGGGSSGL PPSTKKSGQK PVLVDVHAELD RITRKPVTVS PTTPTSPTTEG EAS
Source	Mammalian cell
Target Names	DYNC1LI1
Protein Names	Recommended name: Cytoplasmic dynein 1 light intermediate chain 1 Short name= LIC1 Alternative name(s): Dynein light chain A Short name= DLC-A Dynein light intermediate chain 1, cytosolic
Expression Region	1-523
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.