



Recombinant Human Coiled-coil domain-containing protein 61 (CCDC61)

Product Code	CSB-MP897319HU
Abbreviation	CCDC61
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y6R9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDQPAGLQVD YVFRGVEHAV RVMVSGQVLE LEVEDRMTAD QWRGEFDAGF IEDLTHKTGN FKQFNIFCHM LESALTQSSE SVTLDLLTYT DLESLRNRKM GGRPGSLAPR SAQLNSKRYL ILIYSVEFDR IHYPLPLPYQ GKPDPVVLQG IIRSLKEELG RLQGLDGQNT RDTRENEIWH LREQVSRLAS EKRELEAQLG RSREEALAGR AARQEAEALR GLVRGLELEL RQERGLGHRV AGRRGQDCRR LAKELEEAKA SERSLRARLK TLTSELALYK RGGRGSRGRG RPARPSPSPT GGRALRFDPT AFVKAKERKQ REIQMKQQQR NRLGSGGSGD GPSVSWSRQT QPPAALTGRG DAPNRSRNRS SSVDSFRSRC SSASSCSDLE DFSELSRGG HRRRGKPPSP TPWSGSNMKS PPVERSHHQQ SLANSGGWVP IKEYSSEHQA ADMAEIDARL KALQEYMNRL DMRS
Source	Mammalian cell
Target Names	CCDC61
Protein Names	Recommended name: Coiled-coil domain-containing protein 61
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.