



Recombinant Schizosaccharomyces pombe Uncharacterized WD repeat-containing protein C2A9.03 (SPBC2A9.03)

Product Code	CSB-BP897332SXV
Abbreviation	SPBC2A9.03
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y7K5
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSNGEFIDQY LLPESGAILQ DYLDYFYTS WYTGNRNLIQ LDHYGFPTAS NFELNPKSTD NCGTVDFQGI PWGNKAIGKR KNFRLYRLHT YSLSCNHSDW SPEELSLDTV QVAAEDSFYE FQRFLKNPNC LRNHHHLCKS VKVVGESHSVI YVSNNSIRME NTLHDNSVQI IDLKNNILQK KAISCLDTTN AVGCLGTYG QYCIFNFKNG SFSTRKLSDD SINHIQLFSE KKFILGCNDG AVYVEDVQRT SPILVYRCSN AVNCVQLFDE TPLALVCTDS KDAMLIDIRS NSVSAKLPHS DSVLTGSWHP NGNILATGGQ DTTAKVWDIR ALGKSFRNLG SQMSAVTNLC FSPLGESLAV TEQADVFQLF DTKLFNTSQV IDFFGEIAGT AYSSSGDSL IGVDDPMLGG ILEFRRPLQH PLSFIF
Source	Baculovirus
Target Names	SPBC2A9.03
Protein Names	Recommended name: Uncharacterized WD repeat-containing protein C2A9.03
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.