



Recombinant *Schizosaccharomyces pombe* Uncharacterized WD repeat-containing protein C63.06 (SPCC63.06)

Product Code	CSB-EP897342SXV
Abbreviation	SPCC63.06
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y7T2
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVTASSSVRV SNEDECFIYD IEQLKDNVVV SYSTGSWSCF DKGTLLEIFK VPKAHTNITG IISCDQLNGV ITCGSEGEIH LWDIRSQAKS AVRSWTQQST PFTCIALNKK NQFATGSELT RSLASVQLWD VRSEQKLIRQ WNDAHNDDIT HLQFHPKDNE LLLTGSVDGL VSLLDTTKEE DSTDPEEDPL LHVINHGASI HLAKFVSKKR VMVLHMESEY AMYKLRDKD EKTWSSNELF SIDDLRAELS CSYVINEVST SDKQFCALAF GDFSNHETKF VLVDTSTGEL KKEPTKLER SEEICRAISF DVKNDVYYSG GEDGLLQAFR V
Source	E.coli
Target Names	SPCC63.06
Protein Names	Recommended name: Uncharacterized WD repeat-containing protein C63.06
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.