



Recombinant Schizosaccharomyces pombe Uncharacterized protein C146.02 (SPBC146.02)

| | |
|--------------------------|---|
| Product Code | CSB-EP897350SXV-B |
| Abbreviation | SPBC146.02 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9Y807 |
| Product Type | Recombinant Protein |
| Immunogen Species | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MNGIRKKPAT LQLHNIKTSC FIDGSDDQSD RTRSSSGDST SNSLKLSHHI SDELINAIYK ELKSPVSPWD KNISQHSPLW LEAASNSIGD SSPLLSPISP LKLRFHKYSV SKQLNSVGKQ GFGLSSNHSR YASPSSPARL ARGALGEELS NSQYRLECAA IQKLVAQKRA NRRSSKSTLK NSANDIEFFD IEEVSSRKGL IDRTPSKRNV TSRVSDAYS R LKSAVKGESS LSFKRNRSKS KRKIQLHLGLT EFNGWEHHQC DWL |
| Source | E.coli |
| Target Names | SPBC146.02 |
| Protein Names | Recommended name: Uncharacterized protein C146.02 |
| Expression Region | 1-263 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |