



Recombinant *Xenopus laevis* Paired box protein Pax-5 (pax5)

Product Code	CSB-MP897422XBE
Abbreviation	pax5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9YH95
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MEIHCKHDPF AAMHRHGGVN QLGGVFNAGR PLPDVVRQRI VELAHQGVRP CDISRQLRVS HGCVSKILGR YYETGSIKPG VIGGSKPKVA TPKVVDKIAD YKRQNPTMFA WEIRDRLAE RVCDNDTVPS VSSINRIIRT KVQQPTNQQI PPSNHSIAST GSVTQVSSVS TDSAGSSYSI SGILGITSSN AETNKRKRDE GIQESPIPNG HSLPGRDFLR KQMRGDLFTQ QQLDVLDRVF ERQHYTDIFT STEDIKPEQT EYSAMASLTG GLDEMSSLT SPSSADIGGT VPGPQSYPLV TGRDLASTTL PGYPPHVPPA GQGSYSAPTL TGMVPGSDFS GSPYSHQYS SYNDSWRFPN PGLLGSPYYY SAARGGAPPA TAAAYDRH
Source	Mammalian cell
Target Names	pax5
Protein Names	Recommended name: Paired box protein Pax-5 Short name= XPax-5
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.