



Recombinant Human Hepatocyte nuclear factor 3-beta (FOXA2)

Product Code	CSB-BP897462HU
Abbreviation	FOXA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y261
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLGAVKMEGH EPSDWSSYYA EPEGYSSVSN MNAGLGMNGM NTYMSMSAAA MGSGSGNMSA GSMNMSSYVG AGMSPSLAGM SPGAGAMAGM GGSAGAAGVA GMGPHLSPSL SPLGGQAAGA MGGLAPYANM NSMSPMYGQA GLSRARDPKT YRRSYTHAKP PYSYISLITM AIQQSPNKML TLSEIQWIM DLFPFYRQNN QRWQNSIRHS LSFNDCFLKV PRSPDKPGKG SFWTLHPDSG NMFENGCYLR RQKRFKCEKQ LALKEAAGAA GSGKAAAGA QASQAQLGEA AGPASETPAG TESPSSASP CQEHKRGLLG ELKGTAAAL SPPEPAPSPG QQQQAAAHLG GPPHHPGLPP EAHLKPEHHY AFNHPFSINN LMSSEQHHH SHHHHQPHKM DLKAYEQVMH YPGYGSPMPG SLAMGPVTNK TGLDASPLAA DTSYYQGVYS RPIMNSS
Source	Baculovirus
Target Names	FOXA2
Protein Names	Recommended name: Hepatocyte nuclear factor 3-beta Short name= HNF-3-beta Short name= HNF-3B Alternative name(s): Forkhead box protein A2 Transcription factor 3B Short name= TCF-3B
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific genes such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. This gene has been linked to sporadic cases of maturity-onset diabetes of the young. Transcript variants encoding different isoforms have been identified for this gene.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.