



Recombinant Human Zinc finger protein 30 homolog (ZFP30)

Product Code	CSB-MP897467HU
Abbreviation	ZFP30
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2G7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MARDLVMFRD VAVDFSQEEW ECLNSYQRNL YRDVILENYS NLVSLAGCSI SKPDVITLLE QGKEPWMVVR DEKRRWTLDL ESRYDTKKLF QGKDIYEMNL SQWKVMERIK SCGLEEQESP HEVCFRQVTK TTSEKMPTYR KLTSLPLYQK SHNREKPYEC GECGKAFRVR QQLTFHQRIH TGEKPYECKE CGKAFRQCAH LSRHQRIHTS DKLYECKKCG KIFTCGSDLR VHQRIHIGEK PYECKECKGA FRVRGQLNLH QRIHTGEKPY ECKECKGAFR QYAHLTRHQR LNIAEKCYEC KECGQAFLCS TGLRLHHKLH TGEKPYECKE CGKAFRVRQQ LTLHQRIHTG EKPYDCKECKG KTFSRGYHLT LHQRIHTGEK PYECKECKQKF FRRYSELISH QGIHIGEKPY ECKECKGAFR LFSQLTQHQS IHFGEKPFKC KECEKTFRLL SQLTQHQS IH TGEKPYDCKE CGKAFLHSS LIQHQRHSG EKPYKCKECK KAFRQHSHLT YHQRIHNVT
Source	Mammalian cell
Target Names	ZFP30
Protein Names	Recommended name: Zinc finger protein 30 homolog Short name= Zfp-30 Alternative name(s): Zinc finger protein 745
Expression Region	1-519
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.