



Recombinant Human Krueppel-like factor 13 (KLF13)

Product Code	CSB-EP897476HU-B
Abbreviation	KLF13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2Y9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAAYVDHF AAELVSMSS RAVVHGPREG PESRPEGAAV AATPTLPRVE ERRDGKDSAS LFVVARILAD LNQQAPAPAP AERREGAAAR KARTPCRLPP PAPEPTSPGA EGAAAAPPSP AWSEPEPEAG LEPEREPGPA GSGEPGLRQR VRRGRSRADL ESPQRKHKCH YAGCEKVYVK SSKLKAHLRT HTGERPFACS WQDCNKKFAR SDELARHYRT HTGEKKFSCP ICEKRFMRSD HLTKHARRHA NFHPGMLQRR GGSRTGSL S DYSRSDASSP TISPASSP
Source	E.coli
Target Names	KLF13
Protein Names	Recommended name: Krueppel-like factor 13 Alternative name(s): Basic transcription element-binding protein 3 Short name= BTE-binding protein 3 Novel Sp1-like zinc finger transcription factor 1 RANTES factor of late activated T-lymph
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.