



Recombinant Human Dehydrogenase/reductase SDR family member 7 (DHRS7)

Product Code	CSB-EP897480HU-B
Abbreviation	DHRS7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y394
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	DG DLTLLWAEWQ GRRPEWELTD MVVWVTGASS GIGEELAYQL SKLGVSLVLS ARRVELERV KRRCLENGNL KEKDILVLPL DLTDTGSHEA ATKAVLQEFGR RIDILVNNGG MSQRSLCMDT SLDVYRKLIE LNYLGTVSLT KCVLPHMIER KQGKIVTVNS ILGIISVPLS IGYCASKHAL RGFFNGLRTE LATYPGIIVS NICPGPVQSN IVENSLAGEV TKTIGNNGDQ SHKMTTSRCV RLMLISMAND LKEVWISEQP FLLVTYLWQY MPTWAWWITN KMGKKRIENF KSGVDADSSY FKIFKTKHD
Source	E.coli
Target Names	DHRS7
Protein Names	Recommended name: Dehydrogenase/reductase SDR family member 7 EC=1.1.-.- Alternative name(s): Retinal short-chain dehydrogenase/reductase 4 Short name= retSDR4
Expression Region	29-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.