



Recombinant Homogentisate 1,2-dioxygenase (hgo-1)

Product Code	CSB-MP897510CXY
Abbreviation	hgo-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y041
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSEFDELKYL TGFGNEHATS DPRVPDALPV GQNSPQKCSH GLYAEQLSGT AFTAPRSQNN RSWLYRIRPS VIHRPFAMK ENDQHWTTNF SSIPPNNQY RWNPFPLPTK EGVTFVDNLY TVCGGGDVIS RTGLAIHQFS CNASMEHTAM YNSDGDFLIV PQQGALEITT EFGRLLVNPQ EIAVIPQGIR FSVAVRGPSR GYILEVYGTH FQLPDLGPIG ANGLANPRDF EAPVAWFEDL DVEFTIINKY QGSWFQAKQG HSPFDVVGWH GNYVPYKYDL KKFVINTVS FDHCDPSIFT VLTAPSVKHG TAIADVFIFP PRWGCADNTF RPPYYHRNCM SEYMGLITGC YEAKEGGFKP GGGSLHSMPT PHGPDFNCFE MASNADLPKQ RVAEGTMSFM FESSLNMAIT NWAVYQNVDK DYYKDWQPLK KHFTMPK
Source	Mammalian cell
Target Names	hgo-1
Protein Names	Recommended name: Homogentisate 1,2-dioxygenase EC= 1.13.11.5 Alternative name(s): Homogentisate oxygenase Homogentisic acid oxidase Homogentisicase
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.